

Standard Operation Procedure (S.O.P.) for cRNA Purification using the MEGAclean Kit

S.O.P. prepared by Andrea Rodríguez, 2021

Modified by José E Lizardi, Ph.D., 2022

This S.O.P. will direct you in the procedure of cRNA purification from enzymatic reaction, such as *in vitro* transcription, through the MEGAclean Kit (cat. no. AM1908). This kit can also be employed to purify antisense RNA amplifications (aRNA). The afore mentioned method provides a simple and fast way of efficiently removing nucleotides, short oligonucleotides, proteins, and salts from RNA. The kit is appropriate for single-stranded RNA larger than 100 nt and double-stranded RNA larger than 200 bp, and it recovers from 1 ng to 500 µg of RNA with a 70% or higher efficiency. The purified cRNA could be employed in applications that requires high purity cRNA. In our case, the cRNA will be microinjected in *Xenopus laevis* oocytes.

The Elution Solution's composition is 0.1 mM EDTA (pH 8) in nuclease-free water. However, the eluted cRNA can be directly microinjected in X. laevis oocytes without any other manipulation.

Important:

- **Check if ethanol (100%) was added to the Wash Solution (a box in the side label should be checked). If not, add 20 mL of ethanol before use.**
- **It is recommended to clean the bench and pipettes with an RNase decontamination solution. Use 70% alcohol if decontamination solution is not available.**
- **The Spin Column is made up of two parts: the Filter Cartridge (clear tube with a white silica membrane band) and the Collection Tube (clear).**
- **The Elution Tube is the same than the Collection Tube.**
- **Microcentrifuge must be capable of attaining 10,000–15,000 x g.**
- **Filter Cartridges should not be subjected to over 16,000 x g for risk of mechanical damage and/or glass contamination in the sample.**
- **All tips and microtubes must be autoclaved (RNase free).**
- **All volumes are given on a per culture tube basis.**

Materials:

- 70% Alcohol or RNase Decontaminant (RNaseZap)
- Waste Container with Detergent (1% Alconox, stored at room temperature)
- Micropipettes and Tips for 100, 200, and 1000 µL
- 100% Ethanol
- 70% Ethanol
- 5 M Ammonium Acetate
- Binding Solution (stored at 4 °C)
- Wash Solution (stored at 4 °C)
- Elution Solution (stored at 4 °C)
- Spin Column (Filter Cartridges + Collection Tubes)
- Elution Tube
- Dry Bath at 95 °C
- Microcentrifuge for 1.5 mL
- Vortex

Procedure:

- a) The working area should be sanitized using RNase decontamination solution or 70% ethanol prior to placing any materials on the bench. Prepare the waste container with 1% alconox.
- b) Pre-heat 55 μL of Elution Solution in the dry bath at 95 $^{\circ}\text{C}$.
- c) Add Elution Solution to the cRNA sample to bring the cRNA mixture to 100 μL . Vortex gently.
- d) Add 350 μL of Binding Solution to the cRNA mixture. **Mix gently by pipetting.**
- e) Add 250 μL of 100% Ethanol to the cRNA mixture. **Mix gently by pipetting.**
- f) Transfer the cRNA mixture to the Spin Column and centrifuge at 15,000 $\times g$ for 1 minute.
- g) Remove the Filter Cartridge from the Spin Column and discard the liquid in the Collection Tube into the waste container. Return the Filter Cartridge to the Collection Tube.
- h) Wash the Filter Cartridge by adding 500 μL of Wash Solution and centrifuge the Spin Column at 15,000 $\times g$ for 1 minute.
- i) Discard the Wash Solution as explained above (**step f**) and add 500 μL of Wash Solution to the Filter Cartridge for a second wash.
- j) Centrifuge the Spin Column at 15,000 $\times g$ for 1 minute, discard the Wash Solution as explained above (**step f**), and centrifuge the Spin Column again for 30 seconds at 15,000 $\times g$.
- k) Elute cRNA from the Filter Cartridge as follow,
 - 1) Transfer the Filter Cartridge into a new Elution Tube.
 - 2) Apply 50 μL of **Elution Solution** to the center of the Filter Cartridge. Close the cap of the tube and centrifuge for 1 minute at 15,000 $\times g$ at room temperature to elute the cRNA.
 - 3) *Optional:* Apply another aliquot of 50 μL of **Elution Solution** and repeat the centrifugation to maximize RNA recovery. **Collect the eluate into the same Elution Tube.**
If glass fibers are observed in your sample, they can be removed by centrifuging the sample briefly and then transferring the RNA to a new tube.
- l) **Optional:** Precipitate with 5 M Ammonium Acetate to concentrate cRNA.
 - 1) Add 1:10 volume of 5 M Ammonium Acetate to the purified cRNA (*i.e. If sample was eluted with 50 μL of Elution Solution, 5 μL of 5 M Ammonium Acetate will be added*).
 - 2) Add 100% ethanol at a volume that is 2.5 times that of the elution (*i.e. If sample was eluted with 50 μL of Elution Solution, 137.5 μL of ethanol will be added*), vortex gently, and incubate at -20°C for 30 minutes.
 - 3) Microcentrifuge at 30,000 $\times g$ for 15 min at 4 $^{\circ}\text{C}$ or room temperature, preferably at 4 $^{\circ}\text{C}$.
 - 4) Carefully remove and discard the supernatant.
 - 5) Wash the pellet with 500 μL of cold 70% ethanol, centrifuge again at 30,000 $\times g$, and carefully remove the ethanol.
 - 6) To remove the last traces of ethanol, quickly recentrifuged at 30,000 $\times g$, and aspirate any residual fluid with a very fine tipped pipette or with a syringe needle.
 - 7) Air dry the pellet.
 - 8) Resuspend the pellet using the desired solution and volume
- m) Decontaminate all the equipment that was used with 70% ethanol. Carefully discard the alconox solution in the waste container into the proper waste container and discard the materials in a biological waste (trash can with red bag).

- n) Be sure all equipment are shut down. Return any material and equipment to its proper place.
The pipettes must be returned to their maximal volume.